

p62Dok (Ab-398) Antibody

Cat. #: ANTY021269

Species: Human, Mouse, Rat

Quantity: 100ug

Concentration: 100ug/100ul

Storage and Stability:

Store at -20°C/1 year

Immunogen:

The antiserum was produced against synthesized non-phosphopeptide derived from human p62Dok around the phosphorylation site of tyrosine 398 (E-G-YP-E-L).

Specificity:

p62Dok (Ab-398) antibody detects endogenous levels of total p62Dok protein.

Tested Application: WB and IHC

Application Notes: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Raised In: Rabbit


Purity: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Storage buffer: WB: 1:500~1:1000 IHC: 1:50~1:100

Form: Liquid

Reference:

Michael J. Wick, et, al. (2001) J. Biol. Chem; 276: 42843 - 42850. Paul D. Simoncic, et, al. (2006) Mol. Cell. Biol; 26: 4149 - 4160. Nadia Dubé, et, al. (2004) PNAS; 101: 1834. Adam C. Yopp, et al. (2005) J. Immunol; 175: 2913 - 2924.



For research use only

Antagene, Inc. Toll Free: 1-866-964-2589, Fax: 1-888-225-1868, info@antageneinc.com, www.antageneinc.com