

Monoclonal Antibody to LYN

Cat. #: Mab-607052

Description:

Lyn (also known as p53/56 Lyn) is a membrane-associated protein tyrosine kinase (PTK) mostly expressed in hemopoietic cells which is important in cellular signaling. It contains an SH2 and SH3 domain and has been found to be cleaved after activation of caspases in apoptosis. A member of the Src family of PTKs, there are two known isoforms for Lyn which plays an indispensable role in the Fc epsilon RI (FcεR1) and the B-cell IgM receptor signaling pathway and is essential for Syk activation and Lat phosphorylation after FcεR1 aggregation and can also phosphorylate Tec on multiple residues. Lyn can also be regulated by IL-2 and IL-3. Lyn is a member of the src family of non-receptor protein tyrosine kinases that is predominantly expressed in haematopoietic tissues. Like all members of the src family, lyn is thought to participate in signal transduction from cell surface receptors that lack intrinsic tyrosine kinase activity. It is associated with a number of cell surface receptors including the B cell antigen receptor and immunoglobulin E receptor (FcεRI).

Immunogen/Specificity:

Ni-NTA purified truncated recombinant Lyn expressed in E. Coli strain BL21 (DE3).

Applications :

ELISA: Propose dilution 1: 10,000.

Determining optimal working dilutions by titration test.

Formulation

Antibodies are purified by protein A affinity chromatography

Reference:

1. Sakaguchi, A.Y., et al. Genet. 34: 175.
2. Hibbs, M.L., et al. Biol. 29: 397-400.
3. Williams, J.C., et al. Trends Biochem. Sci. 23: 179-184.

Clone Number: 2H8D7

Isotype: IgG2b

Species: Human

Storage and Stability: stored at -20 C




Figure 1: Western blot analysis using anti-human Lyn monoclonal antibody against truncated Lyn recombinant protein(1), HL-60 cell lysate (2).