



PENINSULA LABORATORIES INC.
DIVISION OF BACHEM
305 OLD COUNTY ROAD, SAN CARLOS, CALIFORNIA 94070
(800) 922-1516 • (650) 592-5392 • FAX (650) 595 4071
www.bachem.com

T-4326 (IHC7100)
Rabbit Anti- Galanin (Human)

Specifications:

- Amount:** Each vial contains 50uL of whole serum in lyophilized form. Reconstitute with 50ul of distilled water for the equivalent of undiluted antiserum.
- Storage:** This polyclonal antibody is supplied as a lyophilized powder. Store at 2-8°C upon receipt. For long-term storage, store at -20°C and can be stable for up to one year. Avoid repeated freeze/thaw cycles.
- Applications:** Immunohistochemistry, Immunocytochemistry, and other applications
- Sample Types:** The antibody can react with a variety of sample preparation including frozen tissue sections, paraffin embedded tissue and fixed cell cultures.

Starting dilution for Immunohistochemical staining:

Immunofluorescence: 1:200

PAP or ABC: 1:500-1:1000

LAB-SA: 1:500

Cross Reaction:

Peptide	%
Galanin (Human)	100
Galanin (Rat)	100
Galanin (Porcine)	100
PHM-27 (Human)	<0.001
Secretin (Human)	0
VIP (Human, Porcine, Rat)	0
Substance P	0

%Cross-Reactivity data is taken from Radioimmunoassay applications