


CATALOGUE #: 8CN4



PRODUCT NAME: Human cystatin c

Source: Pooled human serum.
Material was tested and found negative for HBsAg, HCV, HIV-1 and 2 antibodies, HCV antibodies, and syphilis.

Application: Calibrator or standard in immunoassays

Analysis: Purity > 95 % (SDS-PAGE)

Immunological activity: Immunoreactivity is checked using monoclonal antibodies specific to cystatin C

Purification: Affinity and size-exclusion chromatography

Presentation: Lyophilized from solution in 10 mM Tris-HCl, pH 7.5, 150 mM NaCl
Product is recommended to reconstitute with deionized water

Storage: – 20°C
After reconstitution store at – 20°C

Other information: Avoid repeated freezing and thawing
It is recommended to aliquot the product after reconstitution