

Catalog No. MC2A-TF5

Mouse IgG2a, Isotype control

| | |
|----------------------|---|
| Product type | Isotype control |
| Clone number | HOPC-1 |
| Isotype | IgG2a |
| Buffer | Phosphate buffered saline containing 15ppm ProClin |
| Concentration | 0.5 mg / mL |
| Volume | 100 uL |
| Label | Tide Fluor™ 5WS Ex:649nm, Em:664nm |
| Storage | Store below - 40°C Aliquot to avoid cycles of freeze/thaw. |

RELATED PRODUCTS:

| Product Name | Quantity | Maker | Cat# |
|--|----------|-------|-----------------|
| Mouse IgG2a, Isotype control Label Tide Fluor™ 2WS | 100 uL | CAC | MC2A-TF2 |
| Mouse IgG2a, Isotype control Label Tide Fluor™ 5WS | 100 uL | CAC | MC2A-TF5 |
| Mouse IgG2b, Isotype control Label Tide Fluor™ 2WS | 100 uL | CAC | MC2B-TF2 |
| Mouse IgG2b, Isotype control Label Tide Fluor™ 5WS | 100 uL | CAC | MC2B-TF5 |
| Anti CD9, Human (Mouse) Label Tide Fluor™ 2WS, (Clone :12A12) | 100 uL | CAC | SHI-EXO-M01-TF2 |
| Anti CD9, Human (Mouse) Label Tide Fluor™5WS, (Clone :12A12) | 100 uL | CAC | SHI-EXO-M01-TF5 |
| Anti CD63, Human (Mouse) Label Tide Fluor™ 2WS, (Clone : 8A12) | 100 uL | CAC | SHI-EXO-M02-TF2 |
| Anti CD63, Human (Mouse) Label Tide Fluor™5WS, (Clone : 8A12) | 100 uL | CAC | SHI-EXO-M02-TF5 |
| Anti CD81, Human (Mouse) Label Tide Fluor™ 2WS, (Clone : 12C4) | 100 uL | CAC | SHI-EXO-M03-TF2 |
| Anti CD81, Human (Mouse) Label Tide Fluor™5WS, (Clone : 12C4) | 100 uL | CAC | SHI-EXO-M03-TF5 |

For research use only, Not for diagnostic use.**COSMO BIO Co., LTD.**

【JAPAN】
 TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
 KOTO-KU. TOKYO 135-0016, JAPAN
 Phone: +81-3-5632-9610
 FAX: +81-3-5632-9619
 URL: <https://www.cosmobio.co.jp/>

**COSMO BIO USA**

【Outside Japan】
 2792 Loker Ave West, Suite 101
 Carlsbad, CA 92010, USA
 email: support@cosmobioussa.com
 Phone/FAX: (+1) 760-431-4600
 URL: www.cosmobioussa.com