

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP-AMT-031

Anti AMT

Product type	Primary antibodies
Immunogen	AMT (AAH07546.1, 1 a.a. ~ 289 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MESLVVGDI AELRPNQGTLSLFTNEAGGILDDLIVTNTSEGHLYVVSNAGCWEKDLALMQ DKVRELQNQGRDVGLEVLNALLALQGPTAAQVLQAGVADDLRKLPMFMTSAVMEVFGVSG CRVTRCGYTGEDGVEISVPVAGAVHLATAILKNPEVKLAGLAARDSLRLEAGLCCLYGNDI DEHTTPVEGSLSWTLGKRRRAAMDFPGAKVIVPQLKGRVQRRRVGLMCEGAPMRAHSPIL NMEGTKIGTVTSGCPSPSLKKNVAMGYVPCEYSRPGTMLLVELPSGPCF*
Gene Alias	GCE, GCST, NKH
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~31.90 kDa using AMT transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of AMT expression in transfected 293T cell line by AMT rabbit polyclonal antibody.


Lane 1: AMT transfected lysate(31.90 kDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of AMT transfected lysate.

Immunoprecipitation



Immunoprecipitation of AMT transfected lysate using AMT rabbit polyclonal antibody and Protein A Magnetic Bead, and Immunoblotted with AMT rabbit polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619