

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-CAPN2-011

Anti CAPN2


Product type	Primary antibodies
Immunogen	CAPN2 (AAH21303.1, 1 a.a. ~ 700 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence	MAGIAAKLAKDREAAEGLGSHERAICYLNQDYEARNECLEAGTLFQDPSFPAIPSAALGF KELGPYSSKTRGIEWKRPTEICADPQFIIGGATRTDICQALGDCWLLAAIASLTNEEI LARVVPLNQSFQENYAGIFHFQFWQYGEWVEVVDDRLPTKDGELLFVHSAEGSEFWAL LEKAYAKINGCYEALSGGATTEGFEDFTGGIAEWYELKKPPPNLFKIIQKALQKGSLLGC SIDITSAADSEAITFQKLVKGHAYSVTGAEVEESNGSLQKLIRIRNPWGEVEWTGRWNDN CPSWNTIDPEERERLTRRHEDGEFWMSFSDFLRHYSRLEICNLTPDTLTSDTYKKWKLTK MDGNWRRGSTAGGCRNYPNTFWMNPQYLIKLEEEDEDEEDGESGCTFLVGLIQKHRRRQR KMGEDMHTIGFGIYEVPEELSGQTNHLSKNFFLTNRARERSDTFINLREVLNRFLKPPG EYILVPSTFEPNKDGDGFCIRVFSEKKADYQAVDDEIEANLEEFDISIDDGFRRLFAQ LAGEDAEISAFELQTLRRVLAKRQDIKSDGFSIETCKIMVMDLSDSGSKLGLKEFYIL WTKIQKYQKIYREIDVDRSGTMNSYEMRKALEEAGFKMPCQLHQVIVARFADDQLIDFD NFVRCLVRLETFLFKIFKQLDPENTGTIELDLISWLCFSVL
---------------------------	---

Gene alias	CANP2, CANPL2, CANPml, FLJ39928, mCANP
-------------------	--

Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~80.00 kDa using CAPN2 transfected 293T lysate.
--------------------------	---

Western Blotting (Transfected lysate)



Western Blot analysis of CAPN2 expression in transfected 293T cell line by CAPN2 rabbit polyclonal antibody.

Lane 1: CAPN2 transfected lysate(80.00 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619