

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.


Catalog No. CNP-CASP7-116

Anti CASP7

Product type	Primary antibodies
Immunogen	CASP7 (NP_001218.1, 1 a.a. ~ 303 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MADDQGCIEEQGVEDSANEDSVDAKPDRSSFVPSLFSKSKKKNVTMRSIKTTTRDRVPTYQY NMNFEKLGKCIINNKNFDKVTGMGVRNGTDKDAEALFKCFRSLGFDVIVYNDSCAKMQ DLLKKASEEDHTNAACFACILLSHGEENVIIYGKDGVTPIKDLTAHFRGDRCKTLLEKPKL FFIQACRGTELDGDIQADSGPINDTDANPRYKIPVEADFLFAYSTVPGYYSWRSPGRGSW FVQALCSILEEHGKDLEIMQILTRVNDRVARHFESQSDDPHFHEKKQIPCVVSMLTKELY FSQ
Gene Alias	CMH-1, ICE-LAP3, MCH3
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~34.30 kDa using CASP7 transfected 293T lysate.

IP

Immunoprecipitation



Immunoprecipitation of CASP7 transfected lysate using CASP7 rabbit polyclonal antibody and Protein A Magnetic Bead, and Immunoblotted with CASP7 mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619