

**POLYCLONAL ANTIBODY**

*For research use only. Not for clinical diagnosis.*

**Catalog No. CNP-CCL15-039**

# Anti CCL15

**Product type** Primary antibodies

**Immunogen** CCL15 (NP\_004158.2, 1 a.a. ~ 113 a.a) full-length human protein.

**Form** Liquid  
Protein A purified

**Volume** 100ul

**Specificity** Human

**Host** Rabbit

**Immunogen Sequence**

MKVSVAALSCLMLVAVLGSQAQFINDAETELMMSKLPLENPVVLNSFHFAADCCTSYISQ  
SIPCSLMKSYFETSSECSKPGVIFLTKKGRQVCAPSGPGVQDCMKLKPYSI


**Gene Alias** HCC-2, HMRP-2B, LKN1, Lkn-1, MIP-1d, MIP-5, NCC-3, NCC3, SCY A15, SCYL3, SY15

**Application notes**

**WB**

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~12.20 kDa using CCL15 transfected 293T lysate.

**Western Blotting (Transfected lysate)**



Western Blot analysis of CCL15 expression in transfected 293T cell line by CCL15 rabbit polyclonal antibody.


Lane 1: CCL15 transfected lysate(12.20 KDa).

Lane 2: Non-transfected lysate.

**IP**

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of CCL15 transfected lysate.

**Immunoprecipitation**



Immunoprecipitation of CCL15 transfected lysate using CCL15 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with CCL15 rabbit polyclonal antibody.

**Storage**

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

*For research use only. Not for clinical diagnosis.*



**COSMO BIO Co., LTD.**

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619