

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-COX10-267


Anti COX10

Product type	Primary antibodies
Immunogen	COX10 (AAH00060.1, 1 a.a. ~ 443 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MAASPHLTSSRLLTGCVGGSVWYLERRTIQDSPHKFLHLLRNVNKQWITFQHFSFLKRMV VTQLNRSHNQVRPKPEPVASPFLEKTSSGQAKAEIYEMRPLSPPSLSLRKPNEKELIE LEPDSVIEDSIDVGKETKEEKRWKEMKLQVYDLPGILAQLSKIKLTALVVSTTAAGFALA PGPFDWPCFLTTSVGTGLASCAANSINQFFEVPFDSNMNRTKNRPLVRGQISPLLAVSFA TCCAVPGVAILTLGVNPLTGALGLFNIFLYTCCYTPLKRISIANWVGAVVGAIPPVMGW TAATGSLDAGAFLLGGILYSWQFPHFNALSWGLREDYSRGGYCMMSVTHPGLCRRVALRH CLALLVLSAAAPVLDITTWTFPIMALPINAYISYLGFRFYVDADRRSSRLFFCSLWHLPL LLLLMLTCKRPSGGGDAGPPPS

Gene Alias

Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~48.90 kDa using COX10 transfected 293T lysate.
--------------------------	---

Western Blotting (Transfected lysate)



Western Blot analysis of COX10 expression in transfected 293T cell line by COX10 rabbit polyclonal antibody.

Lane 1: COX10 transfected lysate(48.90 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619