

**POLYCLONAL ANTIBODY**

*For research use only. Not for clinical diagnosis.*

**Catalog No. CNP-DUSP7-079**


# Anti DUSP7

|                           |                                                                                                                                                                                                                                                                                                                                                       |
|---------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product type</b>       | Primary antibodies                                                                                                                                                                                                                                                                                                                                    |
| <b>Immunogen</b>          | DUSP7 (AAH19107.1, 1 a.a. ~ 320 a.a) full-length human protein.                                                                                                                                                                                                                                                                                       |
| <b>Form</b>               | Liquid                                                                                                                                                                                                                                                                                                                                                |
|                           | Protein A purified                                                                                                                                                                                                                                                                                                                                    |
| <b>Volume</b>             | 100ul                                                                                                                                                                                                                                                                                                                                                 |
| <b>Specificity</b>        | Human                                                                                                                                                                                                                                                                                                                                                 |
| <b>Host</b>               | Rabbit                                                                                                                                                                                                                                                                                                                                                |
| <b>Immunogen Sequence</b> | MLRRLRKGNLPIRSIIPNHADKERFATRCKAATVLLYDEATAEWQPEPGAPASVLGLLQ<br>KLRDDGCQAYYLQGGFNKFQTEYSEHCETNVDSSSSPSSSPPTSVLGLGGLRISSDCSDG<br>ESDRELPSSATESDGSVPSSQPAFPVQILPYLYLGCAKDSTNLDVLGKYGIKYILNVTP<br>NLPNAFEHGGFEFTYKQIPISDHWSQNLSQFFPEAISFIDEARSKKCGVLVHCLAGISRSV<br>TVTVAAYLMQKMNLSLNDAYDFVKRKKSNI SPNFNFMGQLLDFERTLGLSSPCDNHASSEQ<br>LYFSTPTNHNLFPLNTLEST |

**Gene Alias** MKP-X, MKPX, PYST2

**Application notes** **WB**  
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~35.31 kDa using DUSP7 transfected 293T lysate.

**Western Blotting (Transfected lysate)**



Western Blot analysis of DUSP7 expression in transfected 293T cell line by DUSP7 rabbit polyclonal antibody.

Lane 1: DUSP7 transfected lysate(35.31 KDa).

Lane 2: Non-transfected lysate.

**Storage**

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

*For research use only. Not for clinical diagnosis.*



**COSMO BIO Co., LTD.**

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619