

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-FARP2-127

Anti FARP2

Product type Primary antibodies

Immunogen FARP2 (AAH21301.1, 1 a.a. ~ 647 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MGEIEGTYRVLQTAGMRLGAQTPVGVTLEPGQTLLPRMQEKHLHLRVKLLDNTMEIFDI
EPKCDGQVLLTQVWKRLNLVECDYFGMEFQNTQSYWIWLEPMKPIIRQIRRPKNVVLRLA
VKFFPPDPGQLQEEYTRYLFALQLKRDLEERLTCADTTAALLTSHLLQSEIGDYDETLD
REHLKVNEYLPQQHCLEKILEFHQKHVGQTPAESDFQVLEIARKLEMYGIRFHMASDRE
GTKIQLAVSHMGVLVFQGTTKINTFNWSKVRKLSFKRKRFLIKLHPEVHGPYQDTLEFLL
GSRDECKNFWKICVEYHTFFRLLDQPKPKAKAVFFSRGSSFRYSGRTQKQLVDYFKDSGM
KRIPYERRHSKTHTSVRALTADLPKQSFPEGLRTPASPSSANAFYSLSPSTLVPSGLP
EFKDSSTSLTDPQVSIVKSPAAERRSGAVAGGPDTPSAQPLGPPALQPGPGLSTKSPQPS
PSSRKSPSLSPAFAQVPLGPAEQGSSPLLSPVLSDAGGAGMDCEEPRHKRVPADAEAYFIV
KEILATERTYLDLEIVTWFRSAVVKEDAMPATLMTLLFSNIDPIYEFHRGFLREVEQR
LALWEGPSKAHTKGSHQRIGDILLRNMRQLKVFQLHEGHVAGVTKME

Gene Alias FIR, FRG, KIAA0793, PLEKHC3

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~73.30 kDa using FARP2 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of FARP2 expression in transfected 293T cell line by FARP2 rabbit polyclonal antibody.

Lane 1: FARP2 transfected lysate(73.30 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619