

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-FOLR1-156


Anti FOLR1

Product type	Primary antibodies
Immunogen	FOLR1 (NP_000793.1, 1 a.a. ~ 257 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MAQRMTTQLLLLLVWVAVVGEAQTRIAWARTELLNVCMNAKHHKEKPGPEDKLHEQCRPW RKNACCSTNTSQAHKDVSILYRFNWNHCGEMAPACKRHFIQDTCLEYECSPNLGPWIIQV DQSWRKERVNLNPLCKEDCEQWWEDCRTSYTCKSNWHKGNWNTSGFNKCAVGAACQPFHF YFPTPTVLCNEIWTHSYKVSNYSRGSGRCIQMWFDPAQGNPNEEVARFYAAAMSGAGPWA AWPFLLSLALMLLWLLS

Gene Alias FBP, FOLR, FR-alpha, MOv18

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~29.80 kDa using FOLR1 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of FOLR1 expression in transfected 293T cell line by FOLR1 rabbit polyclonal antibody

Lane 1: FOLR1 transfected lysate(29.80 KDa).
Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619