

**POLYCLONAL ANTIBODY**

*For research use only. Not for clinical diagnosis.*

**Catalog No. CNP-G0S2-175**


# Anti G0S2

<b>Product type</b>	Primary antibodies
<b>Immunogen</b>	G0S2 (AAH09694.1, 1 a.a. ~ 104 a.a) full-length human protein.
<b>Form</b>	Liquid Protein A purified
<b>Volume</b>	100ul
<b>Specificity</b>	Human
<b>Host</b>	Rabbit
<b>Immunogen Sequence</b>	METVQELIPLAKEMMAQKRKGKMKLYVLGSLALFGVVLGLMETVCSPFTAARRLRDQE AAVAELQAALERQALQKQALQEKGKQQDTVLGGRALSNRQHAS

**Gene Alias** RP1-28O10.2

**Application notes** **WB**  
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~11.44 kDa using G0S2 transfected 293T lysate

**Western Blotting (Transfected lysate)**



Western Blot analysis of G0S2 expression in transfected 293T cell line by G0S2 rabbit polyclonal antibody.


Lane 1: G0S2 transfected lysate(11.44 KDa).

Lane 2: Non-transfected lysate.

## IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of G0S2 transfected lysate.

### Immunoprecipitation



Immunoprecipitation of G0S2 transfected lysate using G0S2 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with G0S2 rabbit polyclonal antibody.

## Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

*For research use only. Not for clinical diagnosis.*



**COSMO BIO Co., LTD.**

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619