



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-GANAB-086

Anti GANAB

Product type	Primary antibodies
Immunogen	GANAB (NP_938148.1, 1 a.a. ~ 944 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence	MAAVAAVAARRRRSWASLVLAFLGVCLGITLAVDRSNFKTCEESSFCKRQRSIRPGLSPY RALDLSLQLGPDSLTVHLIHEVTKVLLVLELQGLQKNMTRFRIDELEPRRPRYRVPDVLV ADPPIARLSVSGRDENSVELTMAEGPYKIILTARPFRLDLEDRLLSVNARGLLFEFH QRAPRVSQGSKDPAEGDGAQPEETPRDGDKPEETQGKAEKDEPGAWREETFKTHSDSKPYG PMSVGLDFSLPGMEHVYGIPEHADNLRKLVTEGGEPYRLYNLDVFQYELYNPMALYGSVP VLLAHNPHRDLGIFWLNAETWVDISSNTAGKTLFGKMMDYLQGSGETPQTDVRWMSETG IIDVFLLLGPSISDVFRQYASLTGTQALPPLFSLGYHQSRWNYRDEADVLEVDQGFDHNDH LPCDVIWLDIEHADGKRYFTWDPSRFPQPRTMLERLASKRRKLVAIVDPHIKVDSGYRVH EELRNLGLYVKTRDGSDEYEGWCWPGSAGYPDFTNPTMRAWWANMFSYDNYEGSAPNLFVW NDMNEPSVFNGPEVTMLKDAQHYGGWEHRDVHNIYGLYVHMATADGLRQRSGGMERPFVL ARAFFAGSQRFGAVWTGDNTAEWDHLKISIPMCLSLGLVGLSFCGADVGGFFKNPEPELL VRWYQMGAYQPFFRAHAHLDTGRREPWLLPSQHNDIIRDALGQRYSLLPFWYTLLYQAGR EGIPVMRPLWVQYPQDVTTFNIDDQYLLGDALLVHPVSDSGAHGVQVYLPQGQGEVWYDIQ SYQKHHGPQTLYLPVTLSIPVFQRRGTIVPRWMRVRSSSECMKDDPITLFFVALSPQGT QGELFLDDGHTFNYQTRQEFLLRRFSFSGNTLVSSSADPEGHFETPIWIERVVIIGAGKP AAVVLQTKGSPESRLSFQHDPETSVLVLRKPGINVASDWSIHLR
--------------------	---


Gene Alias	G2AN, GlulI, KIAA0088
------------	-----------------------

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~106.90 kDa using GANAB transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of GANAB expression in transfected 293T cell line by GANAB rabbit polyclonal antibody.


Lane 1: GANAB transfected lysate(106.90 KDa).

Lane 2: Non-transfected lysate.

IP`

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of GANAB transfected lysate.

Immunoprecipitation



Immunoprecipitation of GANAB transfected lysate using GANAB rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with GANAB purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619