

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-HDAC9-532

Anti HDAC9

Product type Primary antibodies

Immunogen HDAC9 (NP_848512.1, 1 a.a. ~ 1069 a.a) full-length human protein..

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MHSMISSVDVKSEVPVGLPISPLDLRTDLRMMMPVVDPPVREKQLQQEELLIQQQQQIQ
KQLLIAEFQKQHNLTRQHQAQLQEHKLQQELLAIKQQELLEKEQKLEQQRQEQUEVER
HRREQQLPPLRGKDRGRERAVASTEVKQKLQEFLLSKSATKDTPTNGKNHSVSRHPKLWY
TAAHHTSLDQSSPPLSGTSPSYKYTLPGAQDAKDDFPLRKTASEPNLKVRSRLKQKVAER
RSSPLLRRKDGNVVTSFKKRMFEVTESSVSSSSPGSGPSSPNNGPTGSVTENETSVLPPT
PHAEQMVSQQRILIHEDSMNLLSLYTSPSLPNITLGLPAVPSQLNASNSLKEKQKCETQT
LRQGVPLPGQYGGSSIPASSSHPHVTLEGKPPNSSHQALLQHLLLKEQMRQQKLLVAGGVP
LHPQSPLATKERISPGIRGTHKLPRHRPLNRTQSAPLPQSTLAQLVIQQQHQQFLEKQKQ
YQQQIHMNKLLSKSIEQLKQPGSHLEEAEEELQGDQAMQEDRAPSSGNSTRSDSSACVDD
TLGQVGAVKVKEEPVDSDEDAQIQEMESGEQA AFMQPPFLEPTHTRALSVRQAPLAAGVM
DGLEKHRLVSRTHSSPAASVLPHPAMDRPLQPGSATGIAYDPLMLKHQCVCGNSTTHPEH
AGRIQSIWSRLQETGLLNKCEIRIQGRKASLEEIQLVHSEHHSLLYGTNPLDGQKLDPRIL
LGDDSQKFFSSLPCGGLGVDSDTIWNLHSSGAARMAVGCVIELASKVASGELKNGFAVV
RPPGHHAEEESTAMGFCFFNSVAITAKYLRDQLNISKILIVDLVDVHHGNGTQQAFYADPSI
LYISLHRYDEGNFFPGSGAPNEVGTGLGEGYNINIAWTGGLDPPMGDVEYLEAFRTIVKP
VAKEFDPMVLVSAGFDALEGHTPPLGGYKVTAKCFGHLTKQLMTLADGRVVLAEAGHD
LTAICDASEACVNALLGNELEPLAEDILHQSPNMNAVISLQKIIQSKYWKSVRMVAVP
RGCALAGAQLQEETETVSALASLTVDVEQPFAQEDSRTAGEPMEEEPAL


Gene Alias PIG3

**Application
notes**

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~117.60 kDa using HDAC9 transfected 293T lysate...

Western Blotting (Transfected lysate)




Western Blot analysis of HDAC9 expression in transfected 293T cell line by HDAC9 rabbit polyclonal antibody.

Lane 1: HDAC9 transfected lysate(117.60 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of HDAC9 expression in Hela NE (Predicted Mw: 111.2970 kDa).

Western Blotting (Cell lysate)



HDAC9 rabbit polyclonal antibody. Western Blot analysis of HDAC9 expression in Hela NE.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619