

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-IL19-046

Anti IL-19


Product type	Primary antibodies
Immunogen	IL19 (NP_715639.1, 1 a.a. ~ 215 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence
MCTEGAFPHRSACSLPLTHVHTHHVCVPVLWGSVPRGMKLQCVSLWLLGTILILCSVDN
HGLRRCLISTDMHHIEESFQEIKRAIQAKDTFPNVTILSTLETQIIKPLDVCCVTKNLL
AFYVDRVFKDHHQEPNPKILRKISSIANSFLYMQKTLRQCQEQRQCHCRQEATNATRVIH
DYDQLEVHAAAIKSLGELDVFLAWINKNHEVMFSA

Gene Alias IL-10C, MDA1, NG.1, ZMDA1

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~24.60 kDa using IL19 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of IL19 expression in transfected 293T cell line by IL19 rabbit polyclonal antibody.

Lane 1: IL19 transfected lysate(24.60 KDa).
Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619