

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-091


Anti TUBD1

Product type	Primary antibodies
Immunogen	TUBD1 (AAH00258.1, 1 a.a. ~ 398 a.a) full-length human protein..
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MSIVTVQLGQCGNQIGFEVFDALLSDSHSSQGLCSMRENEAYQASCKERFFSEEENGVP ARAVLVDMPEPKVINQMLSKAAQSGQWKYGQHACFCQKQSGSGNNWAYGYSVHGPRHEESIM NIIRKEVEKCDSESGFFIIMSMAGGTGSGLGAFVTQNLEDQYSNSLKMNIQIWPYGTGEV IVQNYNSILTLSHLYRSSDALLHENDAIHKICAKLMNIQISFSDINQVLAHQLGSVFQ PTYSAESSFHYRRNPLGIDRHVWPPLSGLPPLSKMSLNKDLHFNTSIANLVILRGKDVQS ADVEGFKDPALYTSWLKPVNAFNVWKTQRAFSKYEKSAVLVSNSQFLVKPLDMIVGKAWN MFASKAYIHQYTKFGIEEEDFLDSFTSLEQVVASYCNL

Gene Alias FLJ12709, TUBD

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~44.70 kDa using TUBD1 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of TUBD1 expression in transfected 293T cell line by TUBD1 rabbit polyclonal antibody.


Lane 1: TUBD1 transfected lysate(44.70 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of TUBD1 transfected lysate.

Immunoprecipitation



Immunoprecipitation of TUBD1 transfected lysate using TUBD1 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with TUBD1 rabbit polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619