

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-126

Anti KRT34

Product type Primary antibodies

Immunogen KRT34 (AAH41070.1, 1 a.a. ~ 436 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MLYAKPPPTINGIKGLQRKERLKPAHIHLQQLTCFSITCSSTMSYSCCLPSLGCRITSCSS
RPCVPPSCHGYTLPGACNIPANVSNCNWFCESFNGSEKETMQFLNDRLASYLEKVRQLE
RDNAELEKLIQERSQQQEPLLCPSYQSYFKTIEELQQKILCAKAENARLVNIDNAKLAS
DDFRSKYQTEQSLRLLVESDINSIRRIIDELTLCCKSDLESQVESLREELICLKNHEEEV
NTLRSQLGDRNLNVEVDLTAPTVDLNQVLNETRSQYEALVETNRREVEQWFATQTEELNKQV
VSSSEQLQSCQAEIIELRRTVNALEIELQAQHNLRDSENLTESEAHYSSQLSQVQSLI
TNVESQLAEIRCDLERQNQEYQVLLDVRARLECEINTYRSLLSEDCCKLPCNPCATTNAS
GNSCGPCGTSQKGCCN

Gene Alias HA4, Ha-4, KRTHA4, hHa4

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~49.40 kDa using KRT34 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of KRT34 expression in transfected 293T cell line by KRT34 rabbit polyclonal antibody.


Lane 1: KRT34 transfected lysate(49.40 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of KRT34 transfected lysate.

Immunoprecipitation



Immunoprecipitation of KRT34 transfected lysate using KRT34 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with KRTHA4 mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619