

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-135

Anti TUBGCP2

Product type Primary antibodies

Immunogen TUBGCP2 (NP_006650.1, 1 a.a. ~ 902 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MSEFRIHHDVNELLSLLRVHGGDGAEVYIDLLQKNRTPYVTTTSAHSAKVKIAEFSRTP
EDFLKKYDELKSKNTRNLDPLVYLLSKLTEDKETLQYLQQNAKERAELAAAVGSSTTSI
NVPAAASKISMQEELERKQLGSGVATGSTLQQSLELKRKMLRDKQNKKNQSGQHLPIPPAW
VYERPALIGDFLIGAGISTDTALPIGTLPASQESAVVEDLLYVLVGVDGRYVSAQPLAG
RQSRFTFLVDPNLDLSIRELVHRILPVAASYSVTRFIEEKSSFEYGVNHALAAAMRTL
KEHLILVSQLEQLHRQGGLSLQKLWFYIQPAMRTMDILASLATSVDKGECLGGSTLSLLH
DRSFSYTGDSQAQELCLYLTKAASAPYFEVLEKWIYRGIIHDPYSEFMVEEHELKERIQ
EDYNDKYWDQRYTIVQQQIPSFLOKQKADKILSTGKYLNVVRECGHDVTCVPAKEIITLK
ERAYVEQIEKAFNYASKVLLDFLMEEKELVAHLRSIKRYFLMDQGDFVHFMDLAEELR
KPVEDITPPRLEALLELALRMSTANTDPFKDDLKIDLMPHDLITQLLRVLAETKQEKAM
AHADPTTELALSGLEAFSFDYIVKWPLSLIINRKALTRYQMLFRHMFYCKHVERQLCSVWI
SNKTAKQHSLHSAQWFAGFTLRQRMNFVQNIQYYMMFEVMEPTWHILEKNLKSASNID
DVLGHHTGFLDTCLKDCMLTNPELLKVFSKLMSVCMFTNCMQKFTQSMKLDGELGGQTL
EHSTVLGLPAGAEERARKELARKHLAEHADTVQLVSGFEATINKFDKNFSAHLLDLLARL
SIYSTSDCEHGMAVISRLDFNGFYTERLERLSAERSQKATPQVPVLRGPPAPAPRVAVT
AQ


Gene Alias GCP2, MGC138162, SPBC97, Spc97p

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~102.50 kDa using TUBGCP2 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of TUBGCP2 expression in transfected 293T cell line by TUBGCP2 rabbit polyclonal antibody.


Lane 1: TUBGCP2 transfected lysate(102.50 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of TUBGCP2 transfected lysate.

Immunoprecipitation



Immunoprecipitation of TUBGCP2 transfected lysate using TUBGCP2 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with TUBGCP2 rabbit polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619