

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-139

Anti LRSAM1

Product type Primary antibodies

Immunogen LRSAM1 (NP_001005373.1, 1 a.a. ~ 723 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MPLFFRKRKPSEEARKRLEYQMCLAKEAGADDILDISKCELSIIPFGAFATCKVLQKKVL
IVHTNHLTSLLPKSCSLLSLATIKVLDLHDNQLTALPDDLGLTALQVLNVERNQLMQLP
RSIGNLTQLQTLNVKDNKLKELPDTVGE LRSRLTLNISGNEIQRLPQMLAHVRTLEMLSL
DASAMVYPPREVCAGAGTAAILQFLCKESGLEYYPPSQYLLPILEQDGIENS RDSPDGPTD
RFSREELEWQNRFS DYEK RKEQKMLEKLEFERRLELGQREHTQLLQQSSSQKDEILQTVK
EEQSRLEQGLSEHQ RHLNAERQRLQEQLKQTEQNISSRIQKLLQDNQRQKKSSEILKSLE
NERIRMEQLMSITQEETESLRRRDVASAMQQMLTESCKNR LIQMAYESQRQNLVQQACSS
MAEMDERFQQILSWQQMDQNK AISQILQESAMQKAAFEALQVKKDLMHRQIRSIKLIET
ELLQLTQLELKRKSLDTESLQEMISEQRWALSSLLQQLLKEKQQREEELREILTELEAKS
ETRQENYWL IQYQRLLNQKPLSLKLQEEGMERQLVALLEELSAEHYLPIFAHHRLSLDLL
SQMSPGDLAKVGVSEAGLQHEILRRVQELLDAARIQPELKPPMGEVVTPAPQEPPE SVR
PSAPPAELEVQASECVVCLEREAQMIFLNCGHVCCCQCCQPLRTCPLCRQDIAQRLRIY
HSS

Gene Alias

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~83.60 kDa using LRSAM1 transfected 293T lysate


Western Blotting (Cell lysate)



LRSAM1 rabbit polyclonal antibody. Western Blot analysis of LRSAM1 expression in HeLa.

Antibody dilution: 1:1000. Western Blot analysis of LRSAM1 expression in human liver (Predicted Mw: 83.5940 kDa).

Western Blotting (Tissue lysate)




LRSAM1 rabbit polyclonal antibody. Western Blot analysis of LRSAM1 expression in human liver.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of LRSAM1 transfected lysate

Immunoprecipitation



Immunoprecipitation of LRSAM1 transfected lysate using LRSAM1 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with LRSAM1 mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619