

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-141

Anti ARG2

Product type Primary antibodies

Immunogen ARG2 (AAH29050.1, 1 a.a. ~ 354 a.a) full-length human protein.

Form Liquid

Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MSLRGSLSRLLQTRVHSILKKS VHSVAVIGAPFSQGQKRKGVEHGPA AIREAGLMKRLSS
LGCHLKDFGDL SFTPVPKDDLYNNLIVNPRSVGLANQELAEVVSRAVSDGYSCVTLGGDH
SLAIGTISGHARHCPDLCVWVDAHADINTPLTTSSGNLHGQPVSFLLRELQDKVPQLPG
FSWIKPCISSASIVYIGLRDVPPEHFILKNYDIQYFSMRDIDRLGIQKVMERTFDLLIG
KRQRPIHLSFDIDAFDPTLAPATGTPVVGGLTYREGMYIAEEIHNTGLLSALDLVEVNPQ
LATSEEEAKTTANLAVDVIASSFGQTREGGHIVYDQLPTPSSPDESENQARVRI


Gene Alias

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~38.60 kDa using ARG2 transfected 293T lysate.


Western Blotting (Transfected lysate)



Western Blot analysis of ARG2 expression in transfected 293T cell line by ARG2 rabbit polyclonal antibody.

Antibody dilution: 1:1000. Western Blot analysis of ARG2 expression in human placenta (Predicted Mw: 38.5780 kDa).

Western Blotting (Tissue lysate)




ARG2 rabbit polyclonal antibody. Western Blot analysis of ARG2 expression in human placenta.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of ARG2 transfected lysate.

Immunoprecipitation



Immunoprecipitation of ARG2 transfected lysate using ARG2 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with ARG2 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619