



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-147

Anti G3BP

Product type Primary antibodies

Immunogen G3BP1 (NP_005745.1, 1 a.a. ~ 466 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MVMEKPSPLLVGREFVRQYYTLLNQAPDMLHRFYGKNSSYVHGGLD SNGKPADAVYGQKE
IHRKVMSQNFTNCHTKIRHVD AHA TLNDGVVVQVMGLLSNNNQALRRFMQTFVLAPEGSV
ANKFYVHNDIFRYQDEVFGGFVTEPQEESEEEVEEPEERQQTPEVVPDDSGTFYDQAVVS
NDMEEHLEEPVAEPEPDPEPEPEQEPVSEIQEEKPEPVLEETAPEDAQKSSSPAPADIAQ
TVQEDLR TFSWASVTSKNLPPSGAVPVTGIPPHVVKVPASQPRPESKPESQIPPQRPQRD
QRVREQRINIPPQRGPRPIREAGEQGDIEPRRMVRHPDSHQLFIGNLPHEVDKSELKDFF
QSYGNVVELRINSGGKLPNFGFVVFDDSEPVQKVLSNRPIMFRGEVRLNVEEKKTRAARE
GDRRDNR LRGP GPRGGLGGGM RGP RGMVQKPGFGVGRGLAPRQ

Gene Alias G3BP, HDH-VIII, MGC111040

Application notes **IP**
Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of G3BP1 transfected lysate.

Immunoprecipitation



Immunoprecipitation of G3BP1 transfected lysate using G3BP1 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with G3BP monoclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619