

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-161

Anti HMGCR


Product type	Primary antibodies
Immunogen	HMGCR (AAH33692.1, 1 a.a. ~ 835 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence	MLSRLFRMHGLFVASHPWEVIVGTVTLTICMMSNMFTGNNKICGWNYECPKFEEDVLSS DIIILTITRCIAILYIFQFQNLRLQGLSKYILGIAGLFTIFSSFVFSTVVIHFLDKELTG LNEALPFFLLIDLSRASTLAKFALSSNSQDEVRENIARGMAILGPTFTLDALVECLVIG VGTMSGVRQLEIMCCFGCMSVLANYFVFMFFPACVSLVLELSRESREGRPIWQLSHFAR VLEEEENKPNPVTQVRVKMIMSLGLVLVHAHSRWIADPSPQNSTADTSKVSLGLDENVSKR IEPVSLWQFYLSKMISMIDIEQVITLSLALLAVKYIFFEQTETESTLSLKNPITSPVVT QKKVPDNCCRREPMILVRNNQKCDSEEEETGINRERKVEVIKPLVAETDTPNRATFVVGNS SLDTSSVLVTQEPEIELPREPRPNEECLQILGNAEKGAKFLSDAEIQLVNAKHIPAYK LETLMETHERGVSIIRQLLSKKLSEPSLQYLPYRDYNYSLLGGGASSRVLADGMTRGPV VRLPRACDSAIEVKAWLETSEGFAVIKEAFDSTSRFARLQKLHTSIAGRNLIRFQSRSGD AMGMNMISKGTEKALSKLHEYFPEMQILAVSGNYCTDKKPAAINWIEGRGKSVVCEAVIP AKVVREVLKTTTEAMIEVNINKNLVGSAMAGSIGGYNAHAANIVTAIYIACGQDAAQNVG SSNCITLMEASGPTNEDLYISCTMPSIEIGTVGGGTNLLPQQACLQMLGVQGACKDNPGE NARQLARIVCGTVMAGELSLMAALAAGHLVKSHMIHNRSKINLQDLQGACTKKTA
---------------------------	---

Gene Alias

Application notes	IP Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of HMGCR transfected lysate
--------------------------	--

Immunoprecipitation



Immunoprecipitation of HMGCR transfected lysate using HMGCR rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with HMGCR purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619