

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.


Catalog No. CNP- IP-169

Anti BRUNOL4

Product type	Primary antibodies
Immunogen	BRUNOL4 (AAH01946.2, 1 a.a. ~ 484 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MYIKMATLANGQADNASLSTNGLGSSPGSAGHMNGLSHSPGNPSTIPMKDHDAIKLFIGQ IPRNLDEKDLKPLFEEFGKIYELTVLKDRFTGMHKGCAFLTYCERESALKAQSALHEQKT LPGMNRPIQVKPADSESRGGSSCLRQPPSHRKLFGVGMNLKQQSEDDVRRLEAFGNIEEC TILRGPDGNSKGCFAVKYSSHAEQAALHGSQTMPGASSSLVVKFADTDKERTMRRM QQMAGQMGMFNPMAIPFGAYGAYAQALMQQAALMASVAQGGYLNPMAAFAAAQMQQMAA LNMNGLAAAPMTPTSGGSTPPGITAPAVPSIPSPIGVNGFTGLPPQANGQPAAEAVFANG IHPYPAQSPTAADPLQQAYAGVQQYAAAAYPAAYGQISQAFPPPPMIPQQQREGPEGCN LFIYHLPQEFGDAELMQMFLPFGFVSFDNPASAQTAIQAMNGFQIGMKRLKVQLKRPKDA NRPY
Gene Alias	BRUNOL-4, CELF4

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~51.80 kDa using BRUNOL4 transfected 293T lysate

Western Blotting (Transfected lysate)




Western Blot analysis of BRUNOL4 expression in transfected 293T cell line by BRUNOL4 rabbit polyclonal antibody.

Lane 1: BRUNOL4 transfected lysate(51.80 KDa).
Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of BRUNOL4 transfected lysate.

Immunoprecipitation



Immunoprecipitation of BRUNOL4 transfected lysate using BRUNOL4 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with BRUNOL4 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619