

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-177


Anti PDLIM2

Product type	Primary antibodies
Immunogen	PDLIM2 (NP_067643.2, 1 a.a. ~ 352 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MALTVDVAGPAPWGFRTGGGRDFHTPIMVTKVAERGKAKDADLRPGDIIVAINGESAEGM LHAEAQSKIRQSPSPLRLQLDRSQATSPGQTNGDSSLEVLATRFQGSVRTYTESQSSLRS SYSSPTSLSPRAGSPFPSPSSSSSLTGEAAISRSFQSLACSPGLPAADRLSYSGRPGSRQ AGLGRAGDSAVLVLPSPGPRSSRPSMDSEGGSLLLDEDSEVFKMLQENREGRAAPRQSS SFRLQEALEAEERGGTPAFLPSSLSPQSSLPASRALATPPKLHTCEKCSTSIANQAVRI QEGRYRHPGCTCADCGNLKMRGHFWVGDELYCEKHARQRY SAPATLSSRA

Gene Alias FLJ34715, MYSTIQUE, SLIM

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~37.50 kDa using PDLIM2 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of PDLIM2 expression in transfected 293T cell line by PDLIM2 rabbit polyclonal antibody.


Lane 1: PDLIM2 transfected lysate(37.50 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of PDLIM2 transfected lysate.

Immunoprecipitation



Immunoprecipitation of PDLIM2 transfected lysate using PDLIM2 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with PDLIM2 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619