



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-I P-179

Anti ITH1

Product type	Primary antibodies
Immunogen	ITIH1 (NP_002206.1, 1 a.a. ~ 911 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MDGAMGPRGLLLCMYLVSLILQAMPALGSATGRSKSSEKRQAVDTAVDGVFIRSLKVNC KVTSRFAHYVVT SQVVNTANEAREVAFDLEIPKTA FISDFAVTADGNAFIGDIKDKVTAW KQYRKAAISGENAGLV RASGRTMEQFTIHLTVNPQSKVTFQLTYEEVLKRNHMQYEIVIK VKPKQLVHHFEIDVDIFEPQGISKLDAQASFLPKELAAQTIKKSFSGKKGHVLF RPTVSQ QQSCPTCSTSLNGHFKVTYDVSRDKICDLLVANNHFAHFFAPQNL TNMKNVVFVIDIS GSMRGQKV KQTKEALLKILGDMQPGDYFDLVLF GTRVQSWKGS LVQASEANLQAAQDFVR GFSLDEATNLNGG LLRGIEILNQVQESLP ELSNHASILIMLTDGDPTEGVTD RSQILKNV RNAIRGRFPLYNLGFGHNVD FNFLEVMSMENNGRAQRIYEDHDATQQLQG FYSQVAKPLL VDVDLQYPQDAVLALTQNH HKQYYEGSEIVVAGRIADNKQSSFKADVQAHGEGQEFSITC LVDEEEMKKLLRERGHMLENHVERLWAYLTIQELLAKRMKVDREVRANLSSQALRMSLDY GFVTPLTSMSIRGMADQDGLKPTIDKPSEDSPPLEMLGPRRTFVLSALQPSPTHSSSNTQ RLPDRVTGVDTDPHFIIHVPQKEDTLCFNINEEPGVILSLVQDPNTGFSVNGQLIGNKAR SPGQHDGTYFGRLGIANPATDFQLEVTPQNITLNP GFGGPVFSWRDQAVLRQDGVVVTIN KKRNLVSVDDGGT FEVVLHRVWKGSSVHQDFLGFYVLDSHRMSARTHGLLGQFFHPIGF EVSDIHPGSDPTKPDATMVVRNRRLTVTRGLQKDYSKDPWHGAEVSCWFIHNNGAGLIDG AYTDYIVPDIF


Gene Alias H1P, IATIH, ITIH, MGC126415

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~101.40 kDa using ITIH1 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of ITIH1 expression in transfected 293T cell line by ITIH1 rabbit polyclonal antibody.


Lane 1: ITIH1 transfected lysate(101.40 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of ITIH1 transfected lysate..

Immunoprecipitation



Immunoprecipitation of ITIH1 transfected lysate using ITIH1 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with ITIH1 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619