

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-184

Anti NOXA1

Product type	Primary antibodies
Immunogen	NOXA1 (NP_006638.1, 1 a.a. ~ 483 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MASLGDLVRAWHLGAQAVDRGDWARALHLFSGVPAPPARLCFNAGCVHLLAGDPEAALRA FDQAVTKDTCMAVGFFQRGVANFQLARFQEALSDFWLALQLRGHAAIDYTQLGLRFLQ AWEVLHNVASAQCQLGLWTEAASSLREAMSKWPEGSLNGLDSALDQVQRRGSLPPRQVPR GEVFRPHRWHLKHLEPVDLFGKAKVVASAIPDDQGWGVRPQQPQGPANHDARSLIMDSP RAGTHQGPLDAETEVGADRCTSTAYQEQRPPQVEQVGKQAPLSPGLPAMGGPGPGPCEDPA GAGGAGAGGSEPLVTVTVQCAFTVALRARRGADLSSLRALLGQALPHAQLGQLSYLAPG EDGHWWPIPEEESLQRAWQDAAACPRGLQLQCRGAGGRPVLVYQVVAQHSYSAQGPEDLGF RQGDTVLDLCEEPDVPLAVDQAWLEGHCDGRIGIFPKCFVVPAGPRMSGAPGRLPRSQQG DQP
Gene Alias	FLJ25475, MGC131800, NY-CO-31, SDCCAG31, p51NOX
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~51.70 kDa using NOXA1 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of NOXA1 expression in transfected 293T cell line by NOXA1 rabbit polyclonal antibody.


Lane 1: NOXA1 transfected lysate(51.70 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of NOXA1 transfected lysate.

Immunoprecipitation



Immunoprecipitation of NOXA1 transfected lysate using NOXA1 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with NOXA1 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619