



POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-231

Anti BCAT2

Product type	Primary antibodies
Immunogen	ADSL (NP_000017.1, 1 a.a. ~ 484 a.a) full-length human protein..
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MAAAALGQIWARKLLSVPWLLCGPRRYASSSFKAADLQLEMTQKPHKKPGPGGEPLVFGKT FTDHMLMVEWNDKGWGQPRIQPFQNLTLHPASSSLHYSLQLFEGMKAFKGGKDQQVRLFRP WLNMDRMLRSAMRLCLPSFDKLELLECIRRLIEVDKDWVPDAAGTSLYVRPVLIGNEPSL GVSQPTRALLFVILCPVGAYFPGGSVTPVSLLADPAFIRAWVGGVGNYKLGGNYGPTVLV QQEALKRGCEQVLWLYGPDHQLTEVGTMNIFVYWTHEDGVLELVTPLNGVILPGVVRQS LLDMAQTWGEFRVVERTITMKQLLRALEEGRVREVFSGTACQVCPVHRILYKDRNLHIP TMENGPELILRFQKELKEIQYGIRAHWMMFPV
Gene Alias	BCAM, BCT2
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~44.30 kDa using BCAT2 transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of BCAT2 expression in transfected 293T cell line by BCAT2 rabbit polyclonal antibody.

Lane 1: BCAT2 transfected lysate(44.30 kDa).
Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of BCAT2 transfected lysate.

Immunoprecipitation



Immunoprecipitation of BCAT2 transfected lysate using BCAT2 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with BCAT2 rabbit polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619