

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-249

Anti PTS

Product Primary antibodies

type

Immunogen PTS (NP_000308.1, 1 a.a. ~ 145 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit


Immunogen Sequence MSTEGGGRRRCQAQVSRRISSFSASHRLYSKFLSDEENLKLFGKCNNPNGHGHNYKVVVTVH
GEIDPATGMVMNLADLKMYEMEEAIMQPLDHNLDMDVPYFADVSTTENVAVYIWDNLQK
VLPVGVLYKVKVYETDNNIVVYKGE

Gene Alias LAMAN, MANB

Application **WB**

notes Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~16.40 kDa using PTS transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of PTS expression in transfected 293T cell line by PTS rabbit polyclonal antibody.


Lane 1: PTS transfected lysate(16.40 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of PTS transfected lysate.

Immunoprecipitation



Immunoprecipitation of PTS transfected lysate using PTS rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with PTS rabbit polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619