

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-250

Anti SGSH

Product type Primary antibodies

Immunogen SGSH (NP_000190.1, 1 a.a. ~ 502 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MSCPVPACCALLLVGLCRARPRNALLLLADDGGFESGAYNNSAIATPHLDALARRSLLF
RNAFTSVSSCSPSRASLLTGLPQHQNMGMYGLHQDVHHFNSFDKVRSLPLLSQAGVRTGI
IGKKHVGPETVYPFDFAYTEENGSVLQVGRNITRIKLLVRKFLQTQDDRPFFLYVAFHDP
HRCGHSQPQYGTFCFKFGNGESGMGRIPDWTPQAYDPLDVLVPYFVPNTPAARADLAAQY
TTVGRMDQGVLVQLQELRDAGVLNDTLVIFTSDNGIPFPSGRTNLYWPGTAEPPLLVSPE
HPKRWGQVSEAYVSLDLTPTILDWFSIPYPSYAIFGSKTIHLTGRSLLPALEAEPLWAT
VFGSQSHHEVTMSYPMRSVQHRHFRLVHNLNFKMPFPIDQDFYVSPTFQDLLNRTTAGQP
TGWYKDLRHYYYRARWELYDRSRDPHETQNLATDPRFAQLLEMLRDQLAKWQWETHDPWV
CAPDGVLEEKLSPPQCQLHNEL

Gene Alias HSS, MPS3A, SFMD

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~56.70 kDa using SGSH transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of SGSH expression in transfected 293T cell line by SGSH rabbit polyclonal antibody.

Lane 1: SGSH transfected lysate(56.70 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of SGSH expression in human liver (Predicted Mw: 56.6950 kDa).

Western Blotting (Tissue lysate)




SGSH rabbit polyclonal antibody. Western Blot analysis of SGSH expression in human liver.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of SGSH transfected lysate

Immunoprecipitation



Immunoprecipitation of SGSH transfected lysate using SGSH rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with SGSH rabbit polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619