

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-260


Anti NME5

Product type	Primary antibodies
Immunogen	NME5 (NP_003542.1, 1 a.a. ~ 212 a.a) full-length human protein...
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MEISMPPPQIYVEKTLAIKPDIVDKEEEIQDIILRSGFTIVQRRKLRLSPEQCSNFYVE KYGKMFFPNLTAYMSSGPLVAMILARHKAISYWLELLGPNNSLVAKETHPDSLRAIYGTD DLRNALHGSNDFAAAEREIRFMFPEVIVEPIPIGQAAKDYLNLHIMPTLLEGLTELCKQK PADPLIWLADWLLKNNPNKPKLCHHPIVEEYP

Gene Alias NM23-H5, NM23H5

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~24.20 kDa using NME5 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of NME5 expression in transfected 293T cell line by NME5 rabbit polyclonal antibody.


Lane 1: NME5 transfected lysate(24.20 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of NME5 transfected lysate.

Immunoprecipitation



Immunoprecipitation of NME5 transfected lysate using NME5 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with NME5 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619