

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-277

Anti MECR

Product Primary antibodies

type

Immunogen MECR (NP_057095.2, 1 a.a. ~ 373 a.a) full-length human protein..

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit


Immunogen Sequence MWVCSTLWRVRTPARQWRGLLPASGCHGPAASSYSASAEPARVRALVYGHGHPAKVVEL
KNLELAAVRGSDVRVKMLAAPINPSDINMIQGNYGFLPELPAVGGNEGVAQVVAVGSNVT
GLKPGDWVIPANAGLGTWRTEAVFSEEALIQVPSDIPLQSAATLGVPNCTAYRMLMDFEQ
LQPGDSVIQNASNSGVGQAVIQIAAALGLRTINVVRDRPDIQKLSDRKSLGAEHVITEE
ELRRPEMKNFFKDMPQPRALNLCVGGKSSTELLRQLARGGTMTVTYGGMAKQPVVASVSL
IFKDLKLRGFWLSQWKDHPDQFKELILTLCDLIRRGQLTAPACSQVPLQDYQSALEAS
MKPFISSKQILTM

Gene Alias CGI-63, FASN2B, NRB1

Application **WB**

notes Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~40.50 kDa using MECR transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of MECR expression in transfected 293T cell line by MECR rabbit polyclonal antibody.


Lane 1: MECR transfected lysate(40.50 KDa).

Lane 2: Non-transfected lysate.

IP`

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of MECR transfected lysate.

Immunoprecipitation



Immunoprecipitation of MECR transfected lysate using MECR rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with MECR purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619