



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-302


Anti PMVK

Product type	Primary antibodies
Immunogen	PMVK (NP_006547.1, 1 a.a. ~ 192 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MAPLGGAPRLVLLFSGKRKSGKDFVTEALQSRLGADVCAVLRRLSGPLKEQYAEHGLNFQ RLLDTSTYKEAFRKDMIRWGEEKRQADPGFFCRKIVEGISQPIWLVSDTRRVSDIQWFRE AYGAVTQTVRVVALEQSRQQRGWVFTPGVDDAESECLDNFGDFDWWVIENHGVEQRLEEQ LENLIEFIRSRL

Gene Alias HUMPMKI, PMK, PMKA, PMKASE

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~22.00 kDa using PMVK transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of PMVK expression in transfected 293T cell line by PMVK rabbit polyclonal antibody.


Lane 1: PMVK transfected lysate(22.00 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of PMVK transfected lysate

Immunoprecipitation



Immunoprecipitation of PMVK transfected lysate using PMVK rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with PMVK purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619