

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-312

Anti ASMTL

Product type Primary antibodies

Immunogen ASMTL (AAH02508.1, 1 a.a. ~ 621 a.a) full-length human protein..

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MVLCPVIGKLLHKRVVLASASPRRQEILSNAGLRFEVVPSKFKEKLDKASFATPYGYAME
TAKQKALEVANRLYQKDLRAPDVVIGADTIVTVGGLILEKPVDKQDAYRMLSRLSGREHS
VFTGVAIVHCSSKDHQLDTRVSEFYEETKVKFSELSEELLWEYVHSGEPMKAGGYGIQA
LGGMLVESVHGDFLNVVGFPNLHFCKQLVKLYPPRPEDLRRSVKHDSIPAADTFEDLSD
VEGGGSEPTQRDAGSRDEKAEAGEAGQATAEAECHRTRETLPPFPTRLLELIEGFMLSKG
LLTACKLKVFDLLKDEAPQKAADIASKVDASACGMERLLDICAAMGLLEKTEQGYSTET
ANVYLASDGEYSLHGFIMHNNDLTWNLFITYLEFAIREGTNQHHRALGKKAEDLFQDAYYQ
SPETRLRFMRAMHGMTKLTACQVATAFNLSRFSSACDVGCTGALARELAREYPRMQVTV
FDLPDIIELAHFQPPGPQAVQIHFAAGDFFRDPLPSAELYVLCRILHDWPDDKVHKLLS
KVAESCKPGAGLLLVTLLDEEKVAQRALMQSLNMLVQTEGKERSLGEYQCLELHGFH
QVQVHLGGVLDAILATKVAP

Gene Alias ASMTLX, ASMTLY, ASTML

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~68.80 kDa using ASMTL transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of ASMTL expression in transfected 293T cell line by ASMTL rabbit polyclonal antibody.

Lane 1: ASMTL transfected lysate(68.80 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of ASMTL expression in human liver (Predicted Mw: 68.8570 kDa).

Western Blotting (Tissue lysate)




ASMTL rabbit polyclonal antibody. Western Blot analysis of ASMTL expression in human liver.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of ASMTL transfected lysate.

Immunoprecipitation



Immunoprecipitation of ASMTL transfected lysate using ASMTL rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with ASMTL rabbit polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619