

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-342


Anti LIN28

Product type	Primary antibodies
Immunogen	LIN28 (AAH28566.1, 1 a.a. ~ 209 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MGSVSNQQFAGGCAKAAEEAPEEDAARAADPEQLLHGAGICKWFNVRMGFGFLSMTA RAGVALDPPVDVVFVHQSKLHMEGFRSLKEGEAVEFTFKKSAKGLESIRVTGPGGVFCIGS ERRPKGKSMQKRRSKGDRCYNCGGLDHHAKECKLPPQPKKCHFCQSISHMVASCPKLAQQ GPSAQGKPTYFREEEEEIHSPDLLPEAQN

CSDD1, FLJ12457, LIN-28, LIN28A, ZCCHC1

Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~23.10 kDa using LIN28 transfected 293T lysate
--------------------------	--

Western Blotting (Transfected lysate)



Western Blot analysis of LIN28 expression in transfected 293T cell line by LIN28 rabbit polyclonal antibody.


Lane 1: LIN28 transfected lysate(23.10 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of LIN28 transfected lysate.

Immunoprecipitation



Immunoprecipitation of LIN28 transfected lysate using LIN28 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with LIN28 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619