

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-412


Anti MLF1

Product type	Primary antibodies
Immunogen	MLF1 (NP_071888.1, 1 a.a. ~ 268 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MFRMLNSSFEDDPFFSESI LAHRENMRQMIRSFSEPFGRDLLSISDGRGRAHNRRGHNDG EDSLTHTDVSSFQTMDQMVSNMNRNYMQKLERNFGQLSVDPNGHSFCSSSVMTYSKIGDEP PKVFQASTQTRRAPGGIKETRKAMRDSDSGLEKMAIGHHHIDRAHVIKKS KNKKTGDEEV NQEFINMNESDAHAFDEEWQSEVLKYKPGRHN LGNTRMRSVGHENPGSRELKRREKPQQS PAIEHGRRSNVLGDKLHIKGSSSVKSNKK

Gene Alias

Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~30.60 kDa using MLF1 transfected 293T lysate
--------------------------	---

Western Blotting (Transfected lysate)




Western Blot analysis of MLF1 expression in transfected 293T cell line by MLF1 rabbit polyclonal antibody.

Lane 1: MLF1 transfected lysate(30.60 KDa).
Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of MLF1 expression in mouse testis (Predicted Mw: 30.6270 kDa).

Western Blotting (Tissue lysate)




MLF1 rabbit polyclonal antibody. Western Blot analysis of MLF1 expression in mouse testis.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of MLF1 transfected lysate.

Immunoprecipitation



Immunoprecipitation of MLF1 transfected lysate using MLF1 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with MLF1 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619