

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-414

Anti TRAP1

Product type Primary antibodies

Immunogen TRAP1 (NP_057376.1, 1 a.a. ~ 704 a.a) full-length human protein

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MARELRALLWGRRLRPLLAPALAAVPGGKPILCPRRTTAQLGPRRNPAWSLQAGRLFSTQTAEDKEEPLHSIISSTESVQGSTSKHEFQAETKKLLDIVARSLYSEKEVFIRELISNASDALEKLRHKLVSDGQALPEMEIHLQTNAEKGTTITQDTGIGMTQEELVSNLGTIARSGSKAFLDALQNQAEASSKIIGQFGVGFYSAFMVADRVEVYSRSAAPGSLGYQWLSDGSGVFEIAEASGVRTGTKIIHLKSDCKEFSSEARVRDVVTKYSNFVSFPLYLNGRRMNTLQAIWMMDPKDVGEWQHEEFYRYVAQAHDKPRYTLHYKTDAPLNIRSIFYVPDMKPSMFDVSRELGSSVALYSRKVLIQTKATDILPKWLRFRIGVVDSEDIPLNLSRELLQESALIRKL RDVLQQRLIKFFIDQSKKDAEKYAKFFEDYGLFMREGIVTATEQEVKEDIALLRYESSALPSGQLTSLSEYASRMRAGTRNIYYLCAPNRHLAEHSPYYEAMKKKDTEVLFCFEQFDEL TLLHLREFDKKKLISVETDIVVDHYKEEFEDRSPAAECLSEKETEELMAWMRNVLGSRVTNVKVTLRDLTHPAMVTVLEMGAARHFLRMQQQLAKTQEERAQLLQPTLEINPRHALIKKLNQLRAS EPGAQLLV DQIYENAMIAAGLVDDPRAMVGRLNELLVKALERH

Gene Alias HSP75, HSP90L

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~80.00 kDa using TRAP1 transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of TRAP1 expression in transfected 293T cell line by TRAP1 rabbit polyclonal antibody.

Lane 1: TRAP1 transfected lysate(80.00 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of TRAP1 expression in human kidney (Predicted Mw: 80.1100 kDa).

Western Blotting (Tissue lysate)




TRAP1 rabbit polyclonal antibody. Western Blot analysis of TRAP1 expression in human kidney.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of TRAP1 transfected lysate.

Immunoprecipitation



Immunoprecipitation of TRAP1 transfected lysate using TRAP1 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with TRAP1 rabbit polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619