



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-470

Anti RUNX3

Product type Primary antibodies

Immunogen RUNX3 (AAH13362.1, 1 a.a. ~ 429 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MASNSIFDSFPTYSPTFIRD PSTSRRTPPSPAFPCGGGGGKMGENSGALSAQA AVGP GG
RARPEVRSMVDVLADHAGELVRTDSPNFLCSVLPSHWRCNKTL PVA FKVVALGDVPDGT V
VTVMAGNDENYSAELRNASAVMKNQVARFNDLRFVGRSGRGKSFTLTITVFTNPTQVATY
HRAIKVTVDGPREPRRRHRQKLEDQTKPFDPDFGDLERLRMRVTPSTPSPRGSLSTTSHFS
SQPQTPIQGTSELNPFSDPRQFDRSFPTLPTLTESRFPDPRMHYPGAMSAAFPYSATPSG
TSISLSVAGMPATSRFHHTYLPPPYPGAPQNNQSGPFQANPSPYHLYYGTSSGSYQFSMV
AGSSSGGDRSPTRMLASCTSSAASVAAGNLMNPSLGGQSDGVEADGSHSNSPTALSTPGR
MDEAVWRPY


Gene Alias AML2, CBFA3, FLJ34510, MGC16070, PEBP2aC

Application notes

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of RUNX3 transfected lysate

Immunoprecipitation



Immunoprecipitation of RUNX3 transfected lysate using RUNX3 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with RUNX3 mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619