

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-479

Anti CUX1

Product type Primary antibodies

Immunogen CUX1 (AAH25422.1, 1 a.a. ~ 676 a.a) full-length human protein

Form Liquid
Protein A purified

Volume 100ul

Specificity Human


Host Rabbit

Immunogen Sequence MAANVGSMFQYWKRFDLQQLQRELDATATVLANRQDESEQSRKRLIEQSREFKKNTPEDL
RKQVAPLLKSFQGEIDALSKRSKEAAFLNVYKRLIDVPDPVPALDLGQQQLQKLVQRLH
DIETENQKLRETLEEYNKEFAEVKNQEVTIKALKEKIREYEQTLKNQAETIALEKEQKLQ
NDFAEKERKLQETQMSTTSKLEEAHKVQSLQTALEKTRTELFDLKTKYDEETTAKADEI
EMIMTDLERANQRAEVAQREAETLREQLSSANHSLQLASQIQKAPDVAIEVLTRSSLEVE
LAAKEREIAQLVEDVQRLQASLTKLRENSASQISQLEQQLSAKNSTLKQLEEKLGQADY
EEVKKELNILKSMEFAPSEGAGTQDAAKPLEVLLLEKNRSLQSESAAALRISNSDLSGRCA
ELQVRITEAVATATEQRELIARLEQDLSIIQSIQRPDAEGAAEHRLEKIPPIKEATLAF
YGPAAPASGALPEGQVDSLLSISSQRERFRARNQELEAENRLAQHTLQALQSELDLRA
DNVKLFEKIKFLQSYPGRGSGSDDTELRYSSQYEERLDPFSSFSKRERQRKYLSLSPWDK
ATLSMGRVLVLSNMARTIGFFYTLFLHCLVFLVLYKLAWSESMERDCATFCAKKFADHLH
KFHENDNGAAAGDLWQ

Gene Alias CASP, CDP, CDP/Cut, CDP1, COY1, CUTL1, CUX, Clox, Cux/CDP, G
OLIM6, Nbla10317, p100, p110, p200, p75

Application notes **WB**
Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~77.20 kDa using CUX1 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of CUX1 expression in transfected 293T cell line by CUX1 rabbit polyclonal antibody.


Lane 1: CUX1 transfected lysate(77.20 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of CUX1 transfected lysate

Immunoprecipitation



Immunoprecipitation of CUX1 transfected lysate using CUX1 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with CUX1 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619