

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-509

Anti RPA1


Product type	Primary antibodies
Immunogen	RPA1 (NP_002936.1, 1 a.a. ~ 616 a.a) full-length human protein
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MVGQLSEGAIAAIMQKGDNIKPILQVINIRPITGNSPPRYRLMSDGLNTLSSFMLAT QLNPLVEEEQLSSNCVCQIHRFIVNTLKDGRVVILMEVLKSAEAVGVKIGNPVPYNE GLGQPQVAPPAPAASPAASSRPQPQNGSSGMGSTVSKAYGASKTFGKAAGPSLSHTSGGT QSKVVPIASLTPYQSKWTICARVTNKSQIRTWSNSRGEGLFSLELVDESGEIRATAFNE QVDKFFPLIEVNVKYYYFSKGTCLKANKQFTAVKNDYEMTFNNETSVMPCEDDHHLPTVQF DFTGIDDLENKSKDSLVDIIGICKSYEDATKITVRSNNREVAKRNIYLMDSGKVVATL WGEDADKFDGSRQPVLAIKGARVSDFGGRSLSVLSSTIIANPDPEAYKLGRWFDAEGQ ALDGVVISDLKSGGVGGSNTNWKTLYEKSENLGQGDKPDYFSSVATVVYLRKENCMYQA CPTQDCNKKVIDQQNGLYRCEKCDTEFPNFKYRMILSVNIADFQENQWVTCFQESAEAIL GQNAAYLGELKDKNEQAFEEVFQANFRSFIFRVRVKVETYNDSEIRKATVMDVKPVDYR EYGRRLVMSIRRSALM
Gene Alias	HSSB, REPA1, RF-A, RP-A, RPA70

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~68.10 kDa using RPA1 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of RPA1 expression in transfected 293T cell line by RPA1 rabbit polyclonal antibody.


Lane 1: RPA1 transfected lysate(68.10 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of RPA1 transfected lysate..

Immunoprecipitation



Immunoprecipitation of RPA1 transfected lysate using RPA1 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with RPA1 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619