



POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-534

Anti WARS2

Product type Primary antibodies

Immunogen WARS2 (NP_957715.1, 1 a.a. ~ 220 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MALHSMRKARERWSFIRALHKGSAAPALQKDSKKRVFSGIQPTGILHLGNYLGAIESWV
RLQDEYDSVLYSIVDLHSITVPQDPAVLRQSILDMTAVLLACGINPEKSILFQQSQVSEH
TQLSWILSCMVRLPRLQHLHQWKAKTTKQKHDGTVGLLTYPVLQAADILLYKSTHVPVGE
DQVQHMELVQDLAQGFNKKYGEFFVPESILSMCVLVFLT


Gene Alias TrpRS

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~24.80 kDa using WARS2 transfected 293T lysate

Western Blotting (Transfected lysate)



Western Blot analysis of WARS2 expression in transfected 293T cell line by WARS2 rabbit polyclonal antibody.


Lane 1: WARS2 transfected lysate(24.80 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of WARS2 transfected lysate

Immunoprecipitation



Immunoprecipitation of WARS2 transfected lysate using WARS2 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with WARS2 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619