

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- IP-561

Anti POU2F2


Product type	Primary antibodies
Immunogen	POU2F2 (AAH06101.1, 1 a.a. ~ 400 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MVHSSMGAPEIRMSKPLEAEKQGLDSPSEHTDTERNGPDTNHQNPQNKTSPFSVSPTGPS TKIKAEDPSGDSAPAAPLPPQPAQPHLPQAQLMLTGSQLAGDIQQLQLQLVLPVGHHL QPPAQFLLPQAQQSQPGLLPTPNLFQLPQQQTQGALLTSQPRAGLPTQPPKCLEPPSHPEE PSDLEELEQFARTFKQRRIKLGFTQGDVGLAMGKLYGNDFSQTTISRFEALNLSFKNMCK LKPLLEKWLNDAETMSVDSSLPSPNQLSSPSLGFGLPGRRRKKRTSIETNVRFALEKSF LANQKPTSEEILLIAEQLHMEKEVIRVWFCNRRQKEKRINPCSAAPMLPSPGKPASYS PH MVTPQGGAGTLPLSQASSSLSTTAQTPALKAATRLSACQA
Gene Alias	OCT2, OTF2, Oct-2

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~43.50 kDa using POU2F2 transfected 293T lysate

Western Blotting (Transfected lysate)




Western Blot analysis of POU2F2 expression in transfected 293T cell line by POU2F2 rabbit polyclonal antibody.

Lane 1: POU2F2 transfected lysate(43.50 KDa).
Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of POU2F2 expression in mouse kidney (Predicted Mw: 49.4390 kDa).

Western Blotting (Tissue lysate)




POU2F2 rabbit polyclonal antibody. Western Blot analysis of POU2F2 expression in mouse kidney.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of POU2F2 transfected lysate

Immunoprecipitation



Immunoprecipitation of POU2F2 transfected lysate using POU2F2 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with POU2F2 mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C or one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619