

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-MCM10-348

Anti MCM10


Product type	Primary antibodies
Immunogen	MCM10 (NP_060988.3, 1 a.a. ~ 874 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence	MDEEDNLSLLTALLEENESALDCNSEENNFLTRENGEPDAFDELFDADGDGESYTEEAD DGETGETRDEKENLATLFGDMEDLTDEEEVPASQSTENRVLPAPAPRREKTNEELQEELR NLQEQMKALQEQLKVTTIKQTASPARLQKSPEKSPRPPLKERRVQRIQESTCFSALDVP ALPRTKRVARTPKASPPDPKSSSSRMTSAPSQPLQTISRNPSPGITRGQIVGTPGSSGET TQPICVEAFSGLRLRRPRVSSTEMNKKMTGRKLIRLSQIKEKMAREKLEEIDWVTFGVIL KKVTPQSVNSGKTFISIWKLNLDLRDLTQCVSLFLFGEVHKALWKTEQGTVVGILNANPMKP KDGSEEVCLSIDHPQKVLIMGEALDLGTCKAKKKNGEPCTQTVNLRDCEYCQYHVQAQYK KLSAKRADLQSTFSGGRIPKKFARRGTSKERLCQDGFYGGVSSASYAASIAAAVAPKK KIQTTLNLVVGKTNLIQETRQKLGIQKSLSCSEEFKELMDLPTCGARNLKQHLAKAT ASGIMGSPKPAIKSISASALLKQKQRMLEMRRRKSEEIQKRFLQSSSEVESPAVPSSSR QPPAQPPRTGSEFPRLGAPATMTPKLGRGVLEGDDVLFYDESPPPRPKLSALAEAKKLA AITKLRAKGQVLTKNPNISIKKKQKDPQDILEVKERVEKNTMFSSQAEDELEPARKKRRE QLAYLESEEFQILKAKSKHTGILKEAEAEMQERYFEPLVKKEQMEEKMRNIREVKCRVV TCKTCAYTHFKLLETVCVSEQHEYHWHDGVKRFFKCPCGNRSISLDRLPNKHCSNCGLYKW ERDGMLEKKTGPKIGGETLLPRGEEHAKFLNSLK
---------------------------	---

Gene Alias	CNA43, DNA43, MGC126776, PRO2249
-------------------	----------------------------------

Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~98.10 kDa using MCM10 transfected 293T lysate.
--------------------------	---

Western Blotting (Transfected lysate)



Western Blot analysis of MCM10 expression in transfected 293T cell line by MCM10 rabbit polyclonal antibody.

Lane 1: MCM10 transfected lysate(98.10 kDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619