

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-MK8IP2-045

Anti MAPK8IP2

Product type	Primary antibodies
Immunogen	MAPK8IP2 (AAH09940.2, 1 a.a. ~ 443 a.a) full-length human protein.
Form	Liquid
	Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit

Immunogen Sequence

MADRAEMFSLSTFHSLSPPGCRPPQDISLEEFDDDLSEITDDCGLGLSYDS DHCEKDSL
SLGRSEQPHPICSFQDDFQEFEMIDDNEEEDDEEDDEEEEDAEDSAGSPGGRGTGPSAPRD
ASLVYDAVKYTLVVDEHTQLELVSLRRCAGLGH DSEEDSGGEASEEEEAGAALLGGGQVSG
DTSPDSPDLTFSKKFLNVFVNSTSRSSSTESFGLFSCLVNGEEREQTHRAVFRFIPRHPD
ELELDVDDPVLVEAEEDDFWFRGFNMRTGERGVFPAFYAHAVPGPAKDLLGSKRSPCWVE
RFDVQFLGSVEVPCHQNGILCAAMQKIATARKLTVHLRPPASCDLEISLRGVKLSLSGG
GPEFQRCSHFFQMKNISFCGCHPRNSCYFGFITKHPLLSRFACHVVFVSQESMRPVAQSVG
RAFLEYYYQEHLAYACPTEDIYLE


Gene Alias IB2, JIP2, PRKM8IPL

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~49.10 kDa using MAPK8IP2 transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of MAPK8IP2 expression in transfected 293T cell line by MAPK8IP2 rabbit polyclonal antibody.

Lane 1: MAPK8IP2 transfected lysate(49.10 KDa).

Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of MAPK8IP2 expression in HeLa (Predicted Mw: 87.9750 kDa).


Western Blotting (Cell lysate)



MAPK8IP2 rabbit polyclonal antibody. Western Blot analysis of MAPK8IP2 expression in HeLa.

Antibody dilution: 1:1000. Western Blot analysis of MAPK8IP2 expression in human colon (Predicted Mw: 87.9750 kDa).

Western Blotting (Tissue lysate)



MAPK8IP2 rabbit polyclonal antibody. Western Blot analysis of MAPK8IP2 expression in human colon.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619