



POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP-PDXK-252

Anti PDXK

Product type	Primary antibodies
Immunogen	PDXK (NP_003672.1, 1 a.a. ~ 312 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MEEECRVLSIQSHVIRGYVGNRAATFPLQVLGFEIDAVNSVQFSNHTGYAHWKGQVLNSD ELQELYEGLRLNNMNKYDYVLTGYTRDKSFLAMVVDIVQELKQQNPRLVYVCDPVLGDKW DGECSMYVPEDLLPVYKEKVVPLADIITPNQFEALLSGRKIHSQEEALRVMDMLHSMGP DTVVITSSDLPSPQGSNYLIVLGSQRRRNPAHSVVMERIRMDIRKVDVAVFVGTGDLFAAM LLAWTHKHPNNLKVACEKTVSTLHHVLQRTIQCAKAQAGEGVRPSPMQLELRMVQSKRDI EDPEIVVQATVL
Gene Alias	C21orf124, C21orf97, DKFZp566A071, FLJ31940, FLJ37311, MGC15873, MGC31754, MGC52346, PKH, PNK, PRED79
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~35.10 kDa using PDXK transfected 293T lysate.

Western Blotting (Transfected lysate)




Western Blot analysis of PDXK expression in transfected 293T cell line by PDXK rabbit polyclonal antibody.

Lane 1: PDXK transfected lysate(35.10 KDa).
Lane 2: Non-transfected lysate.

Antibody dilution: 1:1000. Western Blot analysis of PDXK expression in A-431 (Predicted Mw: 35.1020 kDa).

Western Blotting (Cell lysate)



PDXK rabbit polyclonal antibody. Western Blot analysis of PDXK expression in A-431.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619