

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-LCMT1-554

Anti LCMT1

Product type Primary antibodies

Immunogen LCMT1 (NP_057393.2, 1 a.a. ~ 334 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MATRQRESSITSCCSTSSCDADDEGVRGTCEDASLCKRFAVSIGYWHDPYIQHFVRLSKE
RKAPENIRGYFARVHGVSQLIKAFRLKTECHCQIVNLGAGMDTTFWRLKDEDLLPSKYFE
VDFPMIVTRKLHSIKCKPPLSSPILEHSEDTLQMDGHILDSKRYAVIGADLRDLSELEE
KLKKCNMNTQLPTLLIAECVLVYMTPEQSANLLKWAANSFERAMFINYEQVNMGDRFGQI
MIENLRRRQCFLAGVETCKSLESQKERLLSNGWETASAVDMMELYNRLPRAEVSRIESLE
FLDEMELLEQLMRHYCLCWATKGGNELGLKEITY


Gene Alias CGI-68, LCMT, PPMT1

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~38.40 kDa using LCMT1 transfected 293T lysate

Western Blotting (Transfected lysate)



Western Blot analysis of LCMT1 expression in transfected 293T cell line by LCMT1 rabbit polyclonal antibody.


Lane 1: LCMT1 transfected lysate(38.40 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of LCMT1 transfected lysate

Immunoprecipitation



Immunoprecipitation of LCMT1 transfected lysate using LCMT1 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with LCMT1 purified mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619