

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP- PTCD3-378

Anti PTCD3


Product type	Primary antibodies
Immunogen	PTCD3 (NP_060422.4, 1 a.a. ~ 689 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MAVVSASVRWLGLRSRLGQPLTGRRAGLCEQARSCRFYSGSATLSKVEGTDVTGIEEVVIP KKKTWDKVAVLQALASTVNRDTTAVPYVFQDDPYLMPASSLESRSFLLAKKSGENVAKFI INSYPKYFQKDIAEPHIPCLMPEYFEPQIKDISEAALKERIELRKVKASVDMFDQLLQAG TTVSLETTNSLLDLLCYGDPSTDYHFQQTGQSEALEEENDETSRRKAGHQFGVTWRA KNNARIFSLMPEKNEHSYCTMIRGMVKHRAEQALNLYTELLNNRLHADVYTFNALIEA TVCAINEKFEEKWSKILELLRHMVAQKVKNLQTFNTILKCLRRFHV FARSPALQVLREM KAIGIEPSLATYHHIIRLFDQPGDPLKRSSFIIYDIMNELMGKRFSPKDPDDDKFFQSAM SICSSLRDLELAYQVHGLLKTGDNWKFIPDQHRNFYYSKFFDLICLMEQIDVTLKWYED LIPSAYFPHSQTMHLLQALDVANRLEVIPKIWKDSKEYGHTFRSDLREEILMLMARDKH PPELQVAFADCAADIKSAYESQPIRQTAQDWPATSLNCIALFLRAGRTQEAWKMLGLFR KHNKIPRSELLNELMDSAKVSNSPSQAEVVELASAFSLPICEGLTQRVMSDFAINQEQQ EALSNTALTSDSDTDSSSDSDSDTSEGK
	DKFZp666K071, FLJ20758

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~78.50 kDa using PTCD3 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of PTCD3 expression in transfected 293T cell line by PTCD3 rabbit polyclonal antibody.

Lane 1: PTCD3 transfected lysate(78.50 kDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619