

POLYCLONAL ANTIBODY

For research use only. Not for clinical diagnosis.

Catalog No. CNP-RAD1-502

Anti RAD1

Product type Primary antibodies

Immunogen RAD1 (NP_002844.1, 1 a.a. ~ 282 a.a) full-length human protein.

Form Liquid
Protein A purified

Volume 100ul

Specificity Human

Host Rabbit

Immunogen Sequence MPLL TQQIQDEDDQYSLVASLDNVRNLSTILKAIHFREHATCFATKNGIKVTVENAKCVQ
ANAFIQAGIFQEFKVQEESVTFRINLTVLLDCLSI FGSSPMPGTLTALRM CYQGYGYPLM
LFLEEGGVTVCKINTQEPEETLDFDFCSTNVINKIILQSEGLREAFSELDMTSEVLQIT
MSPDKPYFRLSTFGNAGSSHLDYPKDSLMEAFHCNQ TQVNRYKISLLKPSTKALVLSCK
VSIRTDNRGFLSLQYMIRNEDGQICFVEYYCCPDEEVPESES


Gene Alias HRAD1, REC1

Application notes

WB

Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~31.80 kDa using RAD1 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of RAD1 expression in transfected 293T cell line by RAD1 rabbit polyclonal antibody.

Lane 1: RAD1 transfected lysate(31.80 KDa).

Lane 2: Non-transfected lysate.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619