

POLYCLONAL ANTIBODY


For research use only. Not for clinical diagnosis.

Catalog No. CNP-RNF25-149

Anti RNF25

Product type	Primary antibodies
Immunogen	RNF25 (NP_071898.2, 1 a.a. ~ 459 a.a) full-length human protein.
Form	Liquid Protein A purified
Volume	100ul
Specificity	Human
Host	Rabbit
Immunogen Sequence	MAASASAAAGEEDWVLPSEVEVLESIYLDELQVIKNGRTSPWEIYITLHPATAEDQDSQ YVCFTLVLPVPAEYPHEVPQISIRNPRGLSDEQIHTILQVLGHVAKAGLGTAMLYELIEK GKEILTDNNIPHGQCVCICLYGFQEKEAFTKTPCYHYFHCHCLARYIQHMEQELKAQGQEQ EQERQHATTQKAVGVQCPVCREPLVYDLASLKAPEPQQPMELYQPSAESLRQQEERKR LYQRQQRGGIIDEAERNRYFISLQQPPAPAEPESAVDVSKGSQPPSTLAAELSTSPAV QSTLPPPLPVATQHICEKIPGTRSNQQRLGETQKAMLDPPKPSRGPWRQPERRHHPKGEC HAPKGRDQTQELPPPEGPLKEPMDLKPEPHSQGVGPPQEKGPGSWQGPPIPRTRDCVRW ERSKGRTPGSSYPRLPRGQAYRPGTRRESLGLESKDGS
Gene Alias	AO7, FLJ13906
Application notes	WB Antibody dilution: 1:1000 ~ 1:2000. Specific band of ~51.20 kDa using RNF25 transfected 293T lysate.

Western Blotting (Transfected lysate)



Western Blot analysis of RNF25 expression in transfected 293T cell line by RNF25 rabbit polyclonal antibody


Lane 1: RNF25 transfected lysate(51.20 KDa).

Lane 2: Non-transfected lysate.

IP

Antibody dilution 1:15 ~ 1:60 for immunoprecipitation of RNF25 transfected lysate.

Immunoprecipitation



Immunoprecipitation of RNF25 transfected lysate using RNF25 rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with RNF25 mouse polyclonal antibody.

Storage

For continuous use, store at 2-8°C for one-two days. For extended storage, store in -20°C freezer. Working dilution samples should be discarded if not used within 12 hours.

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619